

**ATTACHMENT A:**

| Line Item                 | Description  | Qty | Unit | Price |
|---------------------------|--|-----|------|-------|
| 0001                      | Small Steam Electric Sterilizer                            | 1   | EA   |       |
| 0002                      | Install Small Steam Electric Sterilizer                    | 1   | EA   |       |
| 0003                      | Loading Car, Transfer Carriage, & Track Assembly 20x20x38" | 3   | EA   |       |
| 0004                      | Install Loading Cart/Transfer Carriage                     | 3   | EA   |       |
| 0005                      | Install ProConnect software.                               | 3   | EA   |       |
| <b><i>Total Cost:</i></b> |  |     |      |       |

***Literature on the above items must be sent separately from price schedule.***

**LINE ITEM 1:**

- Chamber size must be 20"x20"x38".
- Chamber clearance needs to be 18" for intermediate shelves and 15" for bottom shelf.
- Sterilizer must be capable of running Steam Flush Pressure Pulse (SFPP) cycles.
- Sterilizer door operation must be by foot pedal to allow hands free operation.
- Must have 12 programmable cycles which can be customized.
- Must be made out of solid 316L stainless steel.



**LINE ITEM 3**

- Must have a stationary bottom shelf and adjustable upper selves.
- Assemble loading car must consist of a framework, wheel subassemblies and shelves.
- Must have locking device operating from handle end of transfer carriage.

